

# LS-HKJ ELECTRONIC BALANCE



**DESCRIPTION : High precision digital balance**

Features	
<b>Model</b>	LS-HKJ
<b>Type</b>	Precision Balance
<b>Material</b>	ABS Plastic (Body), Stainless Steel (Pan)
<b>Product Weight</b>	3.75 kg

## Specifications

<b>Display</b>	Single large LCD with backlight, 6 digits					
	Single large LED, 6 digits					
<b>Capacity</b>	1200 g	2200 g	3200 g	5200 g	6200 g	12000 g
<b>Division</b>	0.1 g	0.1 g	0.1 g	0.1 g	0.1 g	0.1 g
<b>Unit</b>	g, ct, oz, dwt, ozt, gsm, tar, t, kg, tm					
<b>Dimension</b>	Body : 580mm (L) x 370mm (W) x 685mm (H)					
	Pan : 160mm x 160mm (max capacity < 3200 g)					
	Pan : 185mm x 165mm (max capacity > 3200 g)					
<b>Power supply</b>	AC/DC adaptor or 6V/1.3Ah rechargeable battery					
<b>Functions</b>	ON / OFF, CAL, PRINT, MODE, ZERO					
	Counting weighing mode					
	Tare full of the range					
<b>Interface</b>	RS232					