

LS-HKT ANALYTICAL BALANCE



DESCRIPTION : High precision digital balance

Features

Model	LS-HKT
Type	Precision Balance
Material	ABS Plastic (Body), Stainless Steel (Pan)
Product Weight	3.0 kg

Specifications

Display	Single large LCD with backlight, 6 digits							
	Single large LED, 6 digits							
Capacity	320 g	620 g	1200 g	2200 g	3200 g	4200 g	5200 g	6200 g
Division	0.01 g	0.01 g	0.01 g	0.01 g	0.01 g	0.01 g	0.01 g	0.01 g
Unit	g, ct, oz, dwt, ozt, gsm, tar, t, kg, tm							
Dimension	Body : 520mm (L) x 340mm (W) x 612mm (H)							
	Pan : Ø 130mm							
Power supply	AC/DC adaptor or 6V/1.3Ah rechargeable battery							
Functions	ON / OFF, CAL, PRINT, MODE, ZERO							
Interface	RS232							